

## Riboprobe Preparation

We use reagents from Promega's Riboprobe System – catalog # depends on polymerase ordered, eg. P1420 (SP6), P1430 (T3), P1440 (T7), P1450 (T3 + T7), P1460 (SP6 + T7))

1. Dry 125 - 250  $\mu\text{Ci}$  (10 - 20  $\mu\text{l}$ )  $^{35}\text{S}$ -UTP\* for each reaction in Speed Vac (use screw top microfuge tubes) (NEG739H Easy Tides,  $^{35}\text{S}$ -UTP from NEN)
2. Set up reaction in tube with dried  $^{35}\text{S}$  as follows:

Nuclease free water	(9.0 - X) $\mu\text{l}$
5X buffer	4.0
DTT	2.0
NTP Mix **	3.0
RNAsin	1.0
DNA	X (approx. 1 $\mu\text{g}$ )
Polymerase (T3, T7 or SP6)	<u>1.0</u>
TOTAL	20.0 $\mu\text{l}$

\*10  $\mu\text{l}$  is minimum amount of  $^{35}\text{S}$ -UTP, and usually gives good labeling. If probe doesn't label well, try 20  $\mu\text{l}$ .

\*\*NTP mix is 2.5 mM each rATP, rGTP and rCTP (1  $\mu\text{l}$  each)

3. Incubate at 37° C for 75 minutes
4. Add 1  $\mu\text{l}$  RQ1 DNase to each tube- incubate 15 min at 37°C
5. Run over Sephadex G50 column equilibrated with TE
6. Collect fractions by monitoring for first radioactive peak with geiger counter.
8. Precipitate probe by adding 1/10th volume 3 M Na acetate and 2.5 volumes EtOH, in - 20°C for 2 hours to overnight. Centrifuge at >12,000g for 10 minutes. Remove supernatant, wash pellet with 70% EtOH, spin 2 minutes, remove supernatant and allow pellet to air dry. Resuspend in appropriate volume (eg. 100  $\mu\text{l}$ ) nuclease free H<sub>2</sub>O with 10mM DTT.
9. Count 1  $\mu\text{l}$  of each collected fraction in scintillation counter.